

A-80 BTE Amplifier Hearing Aid

Specifications:

Max Output 129 ± 4 dB

Peak Gain ≥ 38 dB

Total Harmonic Wave Distortion If $f < 1000$ Hz, $\leq 15\%$
If $f > 1000$ Hz, $\leq 10\%$

Equivalent Input Noise Level ≤ 30 dB

Frequency Range 450Hz~3500Hz

Environment temperature $-40^{\circ}\text{C} \sim 55^{\circ}\text{C}$

Relative Humidity $\leq 95\%$

Atmospheric Pressure 50kPa~106kPa