

## DHS-10 Digital Halogen Moisture Meter, Table Halogen Moisture Meter

### Specifications

The main technical parameter

Capacity : 50 g

Accuracy grade: Second level

The actual scale dividing values: 0.01g

Calibration of scale diving values:10mg

Moisture determination accuracy:  $\pm 0.5\%$  (sample $\geq 10\text{g}$ )

Moisture content determination of readability :0.02%~0.1%(sample $\geq 2\text{g}$ )

Temperature control to adjust tolerance: $\pm 1^{\circ}\text{C}$

Heating temperature range:  $60^{\circ}\text{C}$ ~ $200^{\circ}\text{C}$ (adjusted by  $1^{\circ}\text{C}$ )

Heating time range: 0min~99min(adjusted by 1min)

Analytical method: Manual, Timing, Automatic

Display parameter: Nine kinds

Moisture content determination: 0%~100%

Shell size: 360mm $\times$ 215mm $\times$ 170mm

Net weight: 4~5kg