

DHS-10A Digital halogen Moisture meter, table halogen Moisture meter

Specifications

The main technical parameter

Capacity: 50 g

Accuracy Rate: Second level

Actual scale distribution values: 0.005 g

Calibration of scale diving values: 10mg

Moisture determination accuracy: $\pm 0.5\%$ (sample $\geq 2g$)

Moisture content Determination of readability: 0.02% ~ 0.1% (sample $\geq 2g$)

Temperature control to adjust tolerance: $\pm 1\text{ }^{\circ}\text{C}$

Heating temperature range: (60 ~ 200) $^{\circ}\text{C}$ (modified by 1 $^{\circ}\text{C}$)

Heating time range: 0min ~ 99min (modified by 1min)

Analytical method: Manual, timing, automatic

Display parameter: Nine kinds

Moisture content Determination: 0% ~ 100%

Shell size: 360mm \times 215mm \times 170mm

Net Weight: 4.5 kg