

DHS-16 Digital halogen Moisture meter, table halogen Moisture meter

Specifications:

The main technical parameter

Capacity : 120g

Accuracy grade: Second level

The actual scale dividing values: 0.005g

Calibration of scale diving values:10mg

Moisture determination accuracy: $\pm 0.5\%$ (sample ≥ 2 g)

Moisture content determination of readability :0.02% \square 0.1%(sample ≥ 2 g)

Temperature control to adjust tolerance: $\pm 1^{\circ}\text{C}$

Heating temperature range: $\square 60^{\circ}\text{C}$ $\square 200^{\circ}\text{C}$ (adjusted by 1°C)

Heating time range: 0min \square 99min(adjusted by 1min)

Analytical method: Manual, Timing, Automatic

Display parameter: Nine kinds

Moisture content determination: 0% \square 100%

Shell size: 360mm \times 215mm \times 170mm

Net weight: 4 \square 5kg