

## DHS-16A Digital Halogen Moisture Meter, Table Halogen Moisture Meter

### Specifications

The main technical parameter

Capacity : 120g

Accuracy grade: Second level

The actual scale dividing values: 0.002g

Calibration of scale diving values:10mg

Moisture determination accuracy:  $\pm 0.5\%$  (sample $\geq 2g$ )

Moisture content determination of readability :0.02% $\square$ 0.1%(sample $\geq 2g$ )

Temperature control to adjust tolerance: $\pm 1^{\circ}C$

Heating temperature range:  $\square 60 \square 200 \square ^{\circ}C$ (adjusted by  $1^{\circ}C$ )

Heating time range: 0min $\square$ 99min(adjusted by 1min)

Analytical method: Manual, Timing, Automatic

Display parameter: Nine kinds

Moisture content determination: 0% $\square$ 100%

Shell size: 360mm $\times$ 215mm $\times$ 170mm

Net weight: 4 $\square$ 5kg