

Body temperature range: 32.0°C~42.5°C(89.6°F~108.5°F);
-Accuracy: 32-35.9°C/89.6-96.6°F ±0.3°C(0.5°F)
 36-39°C/96.8-102.2°F ±0.2°C(0.4°F)
 39.1-42.5°C/102.4-108.5°F ±0.3°C(0.5°F)
- Surface temperature range: 0°C~100°C(32°F~212°F);
- Accuracy: ±1.0°C;
-Measurement distance: 5 to 8cm;
- °C/°F exchange;
- Backlight, LCD display;
- Data memory: 32 sets latest record;
-Over range indication: HI or LO for body temp.
-Size:150x75x40mm
-Net weight: 120g(without battery)
-Power: 9V, 1pc*6F22
-Auto data hold;
- Auto power off