

Forehead Infrared Thermometer



Features:

Elegant compact design,
1 second testing,
High temperature alert,
High accuracy,
High quality and reliable performance c/w quality components.

Specifications:

- 1.Body temperature range: 32.0°C~43.0°C;
- 2.Accuracy: body $\pm 0.3^{\circ}\text{C}$ (32.0°C-34.9°C, 42.1°C-43°C); $\pm 0.2^{\circ}\text{C}$ (35.0°C-42°C);
- 3.Measurement: 1 to 3cm for forehead;
- 4.°C/°F exchange;
- 5.Backlight, LCD display;
- 6.Data memory: 32 sets of latest records;
- 7.Over range indication: HI or LO for body temp. ;
- 8.Alarm value setting; beeping alarm;
- 9.Auto power off after 15s;

- 10.Low battery indication;
- 11.Power:2pcs X1.5V AAA battery (not included);
- 12.Thermometer size(L*W*H): 131*36*98mm

