

Digital Thermometer (ECT-5K),home thermometer

This digital thermometer is with fleetness and accurate reading on LCD, and the head looks like cartoon.

Speifications

1. Display range: 32.0~42.0oC/90.0~107.6oF
2. Accuracy: ± 0.1 oC/ ± 0.2 oF
3. Min. scale: 0.1
4. Measurement time (reference only, it differs from people to people):
 - a) 10 ± 1 seconds in rectal
 - b) 17 ± 2 seconds in oral
 - c) 26 ± 3 seconds under arm
5. Beeper function
6. Auto shut-off

7. Flexible tip and waterproof
8. Battery: 1.5V button battery (LR/SR-41)
9. Size: 134 x 22 x 14mm
10. LCD: 20 x7mm
11. N.W.:12.5g
12. Memory: last measuring reading

OEM Services

1. Customer can change the color of the Digital Thermometer
2. Private color box and instruction are available

All Medical Devices Are Manufactured In Conformity To Directive 93/42 EEC.

MEDICAL EQUIPMENT/ MEDICAL DIAGNOSTIC TEST KITS/ HOSPITAL BED/ DIGITAL THERMOMETER/ PILL BOX/
HAMMER/ BLOOD PRESSURE MONITOR/ SPHYGMOMANOMETER/ STETHOSCOPE/ BLOOD PRESSURE METER/
BLOOD PRESSURE/ MEDICAL INSTRUMENT/ BLOOD GLUCOSE METER/ANEROID
SPHYGMOMANOMETER/MERCURY SPHYGMOMANOMETER/BLOOD PRESSURE MONIOR/BPM/BLOOD
PRESSURE METER/ANEROID BP MONITOR/MEDICAL PRODUCT/DIGITAL SPHYGMOMANOMETER/DIGITAL
BLOOD PRESSURE MONIOR

Packing Info

1. size:135x29x18(mm)
2. CTN:360PCS
3. GW:14.5KG
4. GW:17.0KG