


EVM35 3.5" LCD 2X-32X Low Vision Portable Digital Video [Magnifier](#)



Feature :

With the development of living standard ,Weak sight population especially with amblyopia and presbyopia have difficulty in reading during daily life , And traditional magnifier or presbyopic glasses have many defects such as not enough magnification times ,low definition and border distortion...These defects are seriously affect weak sight people's daily life and activity.

Aidnear Portable Digital Eyesightaiding Device is a full color digital&optical magnified device ,It adopts our own hardware /software platform with self owned intellectual property rights ,It provides nearly 10 times magnified image on 3.5" screen with high definition and high contrast,Except the full color mode ,it also provides black and white mode ,You can enjoy the incomparable and joyful ☐

Specification:

Device Dimension: 16*12*8cm

Weight: 0.4kg

LCD size: 3.5 inch TFT LCD 4:3

LCD Definition: QVGA (320x240)

Lightness:250cd/M

Magnification range: 2x-32x

LED Illumination: 3 Level Manual Control

Working time: 3~5 hours

View Mode: 4 models

Battery capacity: 1650~1700mAh

Input voltage: 110V~240V

Charging Time: 3.5~4 hours