

Economic PH TESTER Pocket tester PH-98102

Directions:

The PH-98102 is an ideal bag tester for educational purposes.

Technology parameters > >

Measuring range: 0.00-14.00 PH

Resolution: 0.01 PH

Accuracy: ± 0.1 ph (@ 20 °C), ± 0.2 ph

Operating Temperature: 0 °C ~ 50 °C

Automatic temperature compensation: 0 °C ~ 50 °C

Power supply: 3 × 1.5 V (AG13)

Dimensions: 97 mm × 29 mm × 36 mm

Weight: 50 gram