

## FTP-49 Desk Polaroscope

Polariscope is used to distinguish the homogeneous gem from nonhomogeneous gem and the crystalline aggregate, and test the abnormal birefringence of the gem. Under certain circumstances, use the conoscopic view to view the interference pattern and the uniaxiality.

### Technical parameter

1. polaroid: Diameter is  $\varnothing$ 49mm
2. object stage: sample cell can rotate  $360^\circ$
3. conoscopic view
4. electric power source:  $\varnothing$ 110-240V  $\varnothing$ 50-60HZ  $\varnothing$
5. light source: Super bright LED