

The display device 5 inches high-definition digital display; TV and AV interface display devices
Color Mode Full color, black and white, black and white, yellow on black, black on yellow, yellow with blue characters, blue, yellow letters, white with blue characters, blue white, yellow and blue characters, blue yellow letters, 14 color modes
Magnification 4X-50 times (Connection Television greater magnification)
Freeze Support
Interface Type USB2.0 / AV interfaces
Function keys power button, freeze key, mode key, zoom key, down key
Power supply adapter supply or 1500MAH lithium battery
Dimensions (mm)
Weight (g)

Press the power button 3 seconds to boot off

+ Key to zoom key, stepless zoom, total magnification 20

- Key to zoom out key

M key mode button to switch color effects, full color - & gt; black and white ... - & gt; blue yellow word

Freeze key image freeze, thaw

Features:

1. Look away, look at the past two kinds of viewing mode, distant scenery (such as blackboards text), near the panoramic text

2. Simple operation, a variety of different color effect mode, very convenient low vision, color blindness, the elderly friend to use.

3. With voice prompts, make the operation more simple and convenient

The frozen state can adjust the magnification size, color mode switching more fully clear viewing images

5. Highlight HD digital display, brightness adjustment, easy to use in different environments

6. function keys with different color-coded, and there is sense of touch, accurate operation, simple

7. Connect your TV, large screen size viewing images