

## Forehead Infrared Thermometer



- 1.Body temperature range: 32.0°C~42.9°C (89.6°F~109.2°F);
- 2.Accuracy: body  
± 0.3°C/0.54°F(32.0°C-34.9°C, 42.1°C-42.9°C);  
± 0.2°C/0.4°F(35.0°C-42°C);
- 3.Display resolution: 0.1°C/°F
- 4.Object surface temperature range: 0-100°C;
- 5.Measurement: 3 to 5cm for forehead;
- 6.Unit conversion: switch between °C/°F;
- 7.Backlight, LCD display;
- 8.Data memory: maximum 32 data points;
- 9.High temperature warning;
- 10.Low battery warning;
- 11.Power:2pcsX1.5V AAA battery (not included);