

Body Non-contact Infrared Thermometer HT-820D

Features

1. Set alarm value 38
2. User selectable or
3. Automatic data hold& auto power off
4. Automatic selection range and display resolution 0.1(0.1)
5. Backlight LCD display
6. Precise non-contact measurement
7. Designed for intense & professional use
8. Undoubtedly durable and very reliable

Specifications

1. Human Body/Forehead : 32.0°C~43.0°C (90°F~109.4°F)
2. Accuracy: $\pm 0.3^{\circ}\text{C}$ (0.54°F)
3. Measuring Distance: 5~15cm (2~5.9in)
4. Automatic Stop: 7sec.
5. Resolution: 0.1°C (0.1°F)
6. Response Time: 1sec
7. Operating Temperature: 10°C~40°C (50°F~140°F)
8. Relative Humidity Rate: less than or equal to 85%
9. Power Supply : 1pieces 9V battery