

Measuring Range: Body: 34.0°C~42.9°C(93.2°F-109.2°F)

Surface: 0.0°C~100.0°C(32.0°F~212.0°F)

Accuracy: 0.0°C-33.9°C:±2°C /(32.0°F-93.0°F:±3.6°F)

34.0°C-35.0°C:±0.3°C /(93.2°F-95.0°F:±0.5°F)

35.1°C-39.0°C:±0.2°C /(95.2°F-102.2°F:±0.4°F)

39.1°C-42.9°C:±0.3°C /(102.4°F-109.2°F:±0.5°F)

43.0°C-100.0°C:±2°C /(109.4°F-212.0°F:±3.6°F)

LCD Display, Resolution is 0.1°C(0.1°F)

Beep Alarms and Color Alarms for Slight Fever and High Fever.

°C/°FChangeable

50 Memory Storages

Auto Shut-off: 7 Seconds without Operation

Buzzer on/off

Alarm point be settable

Battery: "AAA" size battery x 2

Battery Life: Appro.5,000 Times

Setting be adjustable for different kinds of people and different operating conditions.

Purple, pink or Blue selectable, Once the order quantity be up to 2000pcs, the color can be changed to what the client want.