

J306U Digital Hearing Aid

Specifications:

Type : Digital and Programmable, 4 channels

Power: Super power BTE

Fitting Range: ≤ 125 dB

- 1) Intelligent Spectrum Detection Processing System
- 2) Wide Dynamic Range Compression (WDRC)
- 3) 4 Channels
- 4) Adaptive Feedback Cancellation
- 5) 2 Memories
- 6) FRONTWAVE™ Directional Microphone
- 7) Adjustment for High & Low Frequency
- 8) Memory Switch Indicator
- 9) Low Battery Indicator
- 10) Maximum Power Output (MPO) Control
- 11) Telecoil Input (optional)