

3 in 1 Folding pH Meter KL0101.digital ph meter with 0-14 range

Key Features

1. Folding design feature to use expediently
2. Replaceable pH electrode
3. Clear LCD display featuring both pH, temperature and humidity
4. The KL-0101 offers increased accuracy when measuring solutions at different temperatures as 5. it automatically measures the temperature and compensates the pH reading accordingly.

Specifications

1. Measuring range
pH: 0.00 ~ 14.00
C: -5 ~ 50C (Environment)
RH: 10% ~ 99%RH (Environment)
2. Resolution: pH: 0.01PH; C: 0.1C; RH: 1%RH
3. Accuracy PH ± 0.1 PH; C: ± 1 C; RH ± 5 %RH
4. Automatic Temperature Compensation: 0°C~50C
5. Power Supply:4×1.5V(AG-13 Button Battery)
6. Operating Temperature: 0C~50C
7. Dimension:156mm×55mm×23mm
8. Weight: 95g