

## **Digital Thermometer (MT501),high-precision thermometer**

This digital thermometer is with fleetness and accurate reading on LCD, and also this product is a soft head hermometer.

### **Specifications**

1. Temperature range: 32.0~42.9°C.
2. Accuracy:  $\pm 0.1^{\circ}\text{C}$  /  $\pm 0.2^{\circ}\text{F}$
3. Least display: 0.1.
4. Measurement time:
  - 1) Mouth measurement:  $60\pm 10\text{s}$
  - 2) Oexter measurement:  $100\pm 20\text{s}$
5. Power off automatically.
6. 1.5V button battery.
7. The product can remember the last time's measurment result.
8. high-precision thermometer

### **Competitive Advantages**

- Good products -- modern designs
- Product Features -- more than 30 different kinds of thermometers
- Product Performance -- longer longevity
- Quality Approvals -- Certificates of CE, FDA, RoHS, ISO9001:13485
- Exellent Service

high-precision thermometer

### **Product Introduction**

1. Sold at home and abroad for many years.
2. Unified using imported raw materials.
3. With the Japanese advanced production equipment.
4. Each production line follow tseting productdure.

high-precision thermometer

### **OEM Services**

1. Customer can change the color of the thermometer
2. Private color box and instruction are available

### **Packing Info**

1. CTN: 62\*32\*35CM
2. 360PCS/CTN

3. NW:12.5 KG