

Digital Thermometer (MT502)

This digital thermometer is with fleetness and accurate reading on LCD.

Specification

1. Temperature range: 32.0~42.9°C.
2. Accuracy: $\pm 0.1^{\circ}\text{C}$ / $\pm 0.2^{\circ}\text{F}$
3. Least display: 0.1.
4. Measurement time:
 - 1) Mouth measurement: $60\pm 10\text{s}$
 - 2) Oexter measurement: $100\pm 20\text{s}$
5. Power off automatically.
6. 1.5V button battery.

Packing Info

1. CTN:62*32*35CM
2. 540PCS/CTN
3. NW:15.0KG