

Medical Infrared Thermometer JXB-178

Specifications:

- 1. Measuring Range: Human Body / Forehead: 32.0 ° C ~ 42.9 C (90 °F ~ 109 F)
- 2. Accuracy: ± 0.3 C (0.54 F)
- 3. Measuring Distance: 5 ~ 8cm (2 ~ 5.9in)
- 4. Automatic Stop: 5sec.
- 5. Resolution: 0.1 C (0.1 F)
- 6. Response Time: 0.5sec
- 7. Operating Temperature: 10 C ~ 40 C (50 F ~ 140 F)
- 8. Relative Humidity Rate: less than or equal to 85%