

Newest High quality BTE digital Hearing aid WK-209

Features

1. Micro BTE and Air Construction design
2. Full Digital circuitry
3. High frequency adjustment
4. Low frequency adjustment
5. Adaptive feedback cancellation
6. Adaptive noise reduction
7. Fixed ear mode
8. Digital volume control
9. Open fit and Normal earhook exchangeable
10. With Memory Function
11. More natural, clear sound and great listening comfort
12. Gradually and comfortably increases sound levels according to your unique needs by digital volume control
13. Durable, long-lasting hearing instrument inside and out
14. Reduces unwanted background noise by adjusting focus using two different channels.

Specification

1. Product size: 60x12x8mm
2. Weight: 3 g
3. Max. saturation sound pressure level(SSPL): 123dB±4dB
4. Full on Acoustic Gain: 39dB
5. Total Harmonic Distortion(THD):
F>1000H ≤15%
F<1000H ≤10%
6. Frequency Response: 334-5026 Hz
7. Equivalent Input Noise: 29 dB
8. Battery: A13
9. Battery Current Drain: 2 mA
10. Estimated Battery Life: 128 Hours
11. Rated Voltage: D.C 1.4V
12. Special Function: 2Channels