

Newest High quality BTE digital Hearing aid WK-302

Features

Discreet thin-wire design
High and low tone control
Left & right ear adjustability
Mini Receiver

Specification

- 1.Product size:43x15x10 mm
- 2.Weight:4 g
- 3.Max. saturation sound pressure level(SSPL):130dB±4dB
- 4.Full on Acoustic Gain:F=1600Hz,45±5dB
- 5.Total Harmonic Distortion(THD):F>1000H ≤15% F<1000H ≤10%
- 6.Frequency Response:450-3500 Hz
- 7.Equivalent Input Noise:≤32 dB
- 8.Battery:AG13/LR44/L1154
- 9.Battery Current Drain:4 mA
- 10.Estimated Battery Life:≥48 Hours
- 11.Rated Voltage:D.C 1.5V