

PH Meter PH-98103

This compact tester features a large easy-to-read LCD. Unlike other pH testers, the pH electrode can be interchanged with any pH electrode of different size and shape. This allows you to use application specific electrodes and saves the cost of a new meter when it is finally necessary to replace an exhausted electrode. The tester is fast, accurate, lightweight. Testing pH has never been easier!

Specifications:

Measuring range \square 0.00-14.00PH

Resolution \square 0.01PH

Accuracy \square ± 0.1 pH(@20°C)

± 0.2 PH

Power Supply \square 3 \times 1.5V(AG13)

Operating Temperature \square 0°C \square 50°C

Dimensions \square 175mm \times 51mm \times 25mm

Weight \square 57g