

Forehead Infrared Thermometer

Features

1. Set alarm if temperature is exceeded
2. Memorization of the 32 last measures
3. Automatic data hold & auto power off
4. Automatic selection range and display resolution 0.1 ° C (0.1 ° F)
5. Backlight LCD display
6. Precise non-contact measurement
7. Designed for intense & professional use
8. Undoubtedly durable and very reliable

Specifications

1. Measuring Range: Human Body / Forehead: 32.0 ° C ~ 42.9 ° C (90 ° F ~ 109 ° F)
2. Accuracy: ± 0.3 ° C (0.54 ° F)
3. Measuring Distance: 5 ~ 8cm (2 ~ 5.9in)
4. Automatic Stop: 5sec.
5. Resolution: 0.1 ° C (0.1 ° F)
6. Response Time: 0.5sec
7. Operating Temperature: 10 ° C ~ 40 ° C (50 ° F ~ 140 ° F)
8. Relative Humidity Rate: less than or equal to 85%
9. Power Supply: DC 3V (2 batteries AA)
10. 5pcs boxes put in a big box, 8pcs big box / carton