

Key Features:

- 1□Adjustable Alarm Setting, notifies the user if the reading result is higher than the set temperature value
- 2□Memorize/records/recall the last 32 measurements
- 3□Measure body temperature regardless of room temperature
- 4□Selectable body temperature and surface temperature measurement
- 5□Precise non-contact measurement
- 6□Designed for intense & professional use
- 7□Selectable °C / °F
- 8□With auto range and resolution of 0.1°C (0.1°F)
- 9□LCD display with backlight
- 10□Undoubtedly durable and very reliable

Specifications:

Measuring Range:

Human Body/Forehead : 32.0°C~42.9°C (90°F~109°F)

Accuracy: ±0.3°C (0.54°F)

Measuring Distance: 5~8cm (2~3.2in)

Automatic Stop: 5sec.

Resolution: 0.1°C (0.1°F)

Response Time: 0.5sec

Operating Temperature: 10°C~40°C (50°F~140°F)

Storage data:Store 32pcs testing data

3 kinds of back light change according different temperature

Relative Humidity Rate: less than or equal to 85%

Power Supply : DC 3V(2 batteries AA)□used more than 10,000times,the life for using whole thermomter is more than 3,000,000times.