

Body Infrared Thermometer with adjustable alarm setting RC003T

Features

1. Adjustable Alarm Setting, notifies the user if the reading result is higher than the set temperature value
2. Memorize/records/recall the last 32 measurements
3. Measure body temperature regardless of room temperature
4. Selectable body temperature and surface temperature measurement
5. Precise non-contact measurement
6. Designed for intense & professional use
7. Selectable °C / °F
8. With auto range and resolution of 0.1°C (0.1°F)
9. LCD display with backlight
10. Undoubtedly durable and very reliable

Specifications

1. Human Body/Forehead : 32.0°C~42.9°C (90°F~109°F)
2. Accuracy: ±0.3°C (0.54°F)
3. Measuring Distance: 5~8cm (2~3.2in)
4. Automatic Stop: 5sec.
5. Resolution: 0.1°C (0.1°F)
6. Response Time: 0.5sec
7. Operating Temperature: 10°C~40°C (50°F~140°F)
8. Storage data:Store 32pcs testing data
9. 3 kinds of back light change according different temperature
10. 3 kinds of language function(English/French/Spanish)
11. Relative Humidity Rate: less than or equal to 85%
12. Power Supply : DC 3V(2 batteries AA),used more than 10,000times,the life for using whole thermomter is more than 3,000,000times.

Packing Info

1. Unit size:160x82x42mm
2. Box size:182(L)x132(W)x70(H)mm
3. 30pcs/ctn,11.3kgs