

RC004T IR talking thermometer,medical infrared thermometer,body infrared thermometer,infrared forehead thermometer,medical thermometer



Key Features:

Pet Thermometer Electronic

Precision:Measuring tolerance is no more than $\pm 0.30^{\circ}\text{C}$.

Swiftness:Measuring time is less than 1 second.

Easy operation:One key measurement, operate conveniently.

Non-contact:Forehead measurement, no skin contact.

Large screen display:Large LCD screen display, white back light, makes clear digital performance in any ray.

High temperature alarm: Free in temperature alarm setting.

Others:High temperature alarm: Free in temperature alarm setting.Voice function,As much as 32 groups of data storage

Specifications:

1. Normal conditions of use

● Operating temperature: $10^{\circ}\text{C} \sim 40^{\circ}\text{C}$ ($50^{\circ}\text{F} \sim 104^{\circ}\text{F}$)

● Humidity rate: $\leq 85\%$

2. Power: DC 3V (2 batteries AAA)

3. Size: 141 x 44 x 52 mm - 5.6x1.4x2.1in (L x W x H)

4. Weight: 110g

5. Display Resolution: 0.1°C (0.1°F)

6. Measuring range:

● In body mode: $32^{\circ}\text{C} - 42.9^{\circ}\text{C}$ ($89.6^{\circ}\text{F} - 109.2^{\circ}\text{F}$)

● In Surface Temp mode: $0^{\circ}\text{C} - 60^{\circ}\text{C}$ ($32^{\circ}\text{F} - 140^{\circ}\text{F}$)

7. Precision: From 36°C to 39°C (96.8°F - 102.2°F) = $\pm 0.2^{\circ}\text{C}$ \square 0.4°F \square

8. Consumption: $\leq 300\text{mW}$

9. Accuracy: $\pm 0.3^{\circ}\text{C}$ (0.6°F)

10. Measuring distance: 3-5cm(1.2in-2in)

11. Automatic stop: 5 sec.



Payment Methods:

We accept TT bank wire transfer, Western Union, Paypal and escrow payment offered by Alipay Escrow and Escrow.com.

For latter three payment methods, extra charges shall be applied.

Shipping Methods:

To be shipped by express mail, by air or by sea based on customer's request and weight / volume of goods.

Warranty:

We provide a limited quality warranty of 12 months for the items sold dating from their delivery date.