

RC008 CE Approval non-contact Infrared Thermometer

Specifications:

- 1□Adjustable Alarm Setting, notifies the user if the reading result is higher than the set temperature value
 - 2□Memorize/records/recall the last 32 measurements
 - 3□Measure body temperature regardless of room temperature
 - 4□Selectable body temperature
 - 5□Precise non-contact measurement
 - 6□Designed for intense & professional use
 - 7□Selectable °C / °F
 - 8□With auto range and resolution of 0.1°C (0.1°F)
 - 9□LCD display with backlight
 - 10□Undoubtedly durable and very reliable
 - 11)Can also be used to measure temperature of things such as baby bottle, bath/liquid, etc. and area such as room temperature, etc.
- CE Approval non-contact Infrared Thermometer is your best choice.