

RS13A CE & amp; FDA Approval 2013 newest Hearing aid

1. Latest CLASS-D integrated circuit
2. Ultra- low equivalent input noise (EIN)
3. Rich high frequency to make sound clearer
4. High / low frequency control
5. Telecoil input
6. Low battery consumption
7. MAX-SSPL90: 128dB SPL
8. HFA-SSPL90: 120dB SPL
9. Peak Gain: 54dB
10. HFA-FOG: 48dB
11. Frequency Range: 200-5300Hz
12. THD at: 500Hz: 3.8%
800Hz: 1.2%
1600Hz: 0.2%
13. Eq.Input Noise Level: 28dB
14. Battery Current Drain: 0.97mA
15. HFA-SPLITS: 93dB