

RZ8808A Non-contact Body Thermometer

Specifications

- * Quick measurement: 0.5 second
- * Accurate: 0.2 C (0.4F)
- * Non-contact measurement
- * Body mode: 34 to 42.9C ; Surface mode: 0 to 100 C
- * Three color backlight indication: Green, Yellow, Red
- * Green: Safe, Yellow: Slight fever, Red: High fever
- * C / F switchable
- * 50 memories
- * Auto power off 18 seconds
- * Power supply: 2* AA 1.5V batteries
- * Multiple functions, not only for body temperature, but also for tempertaure of room, Object, milk, food, bath water etc.