

Features:

With the improvement of living standards, .weak pseudo-vision people in particular have amblyopia and hyperopia difficult to read in the daily life and traditional magnifying glass or reading glasses have many flaws, such as not enough magnification, low resolution and frame deformation. These defects seriously affect the visually impaired people's daily life and activities.

Aidnear Eyesightaiding portable digital device is a full-color digital amplifier and; optical amplification equipment, using our own hardware / software platform, with its own intellectual property, which provides nearly 10 times larger image on 3.5 "screen, high-definition , high contrast, except for the full-color mode, it also provides a black and white mode, you can enjoy the unparalleled pleasure.

Specifications:

Product Dimensions: 153 × 80 mm

Weight: 0.25 kg

LCD definitions: QVGA (800X480)

Brightness: To 250cd / M

Magnification range: 4X to 32X

Working Hours: 3 ~ 5 hours

Browse mode: four models

Battery capacity: 2500mAh high capacity rechargeable lithium battery

Input voltage: 110V~240V

Charging time: 3.5~four hours