

## YS008 Portable Video Aids Video Magnifier

### Product Features:

1. Fast algorithm design, smooth reading
2. New design, the use of human gravity feel
3. Small size zero distance design, support handwriting, easy to carry
4. Replace the traditional reading glasses
5. To solve low vision groups
6. Long working hours
7. High-end wide-angle lens

### Product parameters:

Magnification Range 2x-25x

LCD screen 3.5-inch high-definition

Display mode Full color / black on white / black on white

Image Sensor 1/6.5" 640X480

Lens 64 degrees

Image freeze Freeze / Release

Video output 50/60HZ

External TF card Maximum support 32G

External TV standard PAL/NTSC

Automatic shut-down no further operation of button up to 5 min, it turns off automatically

Product weight 130g

Battery AAAX3 lithium battery

Continuous use of time 2.5 hours continual use

Accessories video cable, USB cable, carrying bag, Cleaner Wipe Cloth, Instructions

Product Dimensions 128x77x36mm