

Video Microscope M35

Specifications

Image Sensor: 5 Mega Pixels

Lens : extremely clear microscope lens 5G + IR CUT

Focus Range: Manual focus from 10mm to infinity

Filming resolution : 720P / 1080P, with recording function

Photo image resolution : 12 million pixels, 10 million pixels, 8,000,000 pixels, 5,000,000 pixels, 3,000,000 pixels, 2 million-pixel high-definition HD pixels, 1.3 million pixels, 30 million pixels.

Light source: 8 LED light

TF card storage function : up to 32GB

Download Support : USB port (USB connection to download files)

Support Computer Interface: USB2.0 / USB3.0