

## **Newest High quality BTE Hearing aid WK-022S**

### **Features**

Ergonomic design  
Simple and effective to use  
Power ON/OFF switch  
Left & right ear adjustability  
Detachable Receiver

### **Specification**

- 1.Product size:43x17.5x10 mm
- 2.Weight:7 g
- 3.Max. saturation sound pressure level(SSPL):130dB±4dB
- 4.Full on Acoustic Gain:F=1600Hz,45±5dB
- 5.Total Harmonic Distortion(THD):F>1000H ≤15% F<1000H ≤10%
- 6.Frequency Response:450-3500 Hz
- 7.Equivalent Input Noise:≤32 dB
- 8.Battery:AG13/LR44/L1154
- 9.Battery Current Drain:4 mA