

WK-191 BTE Hearing aid

Specifications

1. MaxSSPL ± 4 dB 130
2. Full-on-gain ± 5 dB 46
3. Equivalent input noise level \pm dB ≤ 28
4. Harmonic distortion % ≤ 4
5. Power supply electric voltage \pm V 1.2
6. The sum settles electric current \pm mA 1.6-2.5
7. Continuoue wrok time(H) ≥ 10
8. Battery Capacity(aAH) 40
9. accessories: Earplug:6pcs Hairbrush:1pcs
10. Charger:1pcs Stand charge:1pcs