

1. Not for otitis media or ear maturation.

2. Should be consulting a qualified audiology or otology doctor as following case:

1) Obvious congenital or traumatic ear canal malformation

2) It has active ear discharging or the secretion of history during last 90 days

3) It has suddenly or progressive hearing loss during last 90 days

4) Acute or chronic vertigo

5) Air-bone gap is equal to or greater than 15dB at pure tone 500Hz, 1000Hz and 2000Hz

6) Cerumen embolization or ear canal with foreign body.

7) Ear pain or discomfort

Technical

size: 53*26.5*8.6mm

weight: 5g

working voltage: DC1.5

Working Current: ≤ 1.5 mA

Frequency Response: 500-3200Hz

Maximum Saturation Acoustic Output: 129dB ± 3 dB

Full-on Acoustic Gain: 36dB ± 5 dB

Total Harmonic Distortion : $< 10\%$

Equivalent Input Noise: < 32 dB