

Specification:

Max. OSPL90 (dBSPL) : 128 (± 4 dB)
Max. Full-on Gain (dB) : 40 (± 5 dB)
Frequency Range (Hz) : 180-4200 (± 6 dB)
Total Harmonic Distortion (THD): 800Hz 6% ($\leq 10\%$)
Equivalent Input Noise (dBSPL) : ≤ 26
Battery Current (mA) : 1.6 (± 0.3 mA)
Working Voltage (V) : 1.2
Charging Time (H) : 4-8
Working Time (H) : 15-18

Payment Methods:

We accept TT bank wire transfer, Western Union, Paypal and escrow payment offered by Alipay Escrow and Escrow.com.

For latter three payment methods, extra charges shall be applied.

Shipping Methods:

To be shipped by express mail, by air or by sea based on customer's request and weight/volume of goods.

Warranty:

We provide a limited quality warranty of 12 months for the items sold dating from their delivery date.

[Meter8](#) is your best choice !